


# PULPDENT®

Live in comfort • Smile with confidence



Calcium Hydroxide Paste placed directly into the canal using 27-gauge x 25mm Endo Irrigation Needles without clogging

# TempCanal™ enhanced

## Calcium Hydroxide Canal Treatment Paste


*For routine use between visits and for extended treatment of complicated cases*

## Easier Placement Better Clinical Results

Flows through 27-gauge irrigation needles with 2-side-vents for controlled placement and uniform coating of canal walls, and with closed end to prevent overfills.

- Easy flow formula – will not clog 27-gauge needles
- Non-drying paste – long working time
- pH>12 – proven benefits of high pH material
- Non-setting – easily removed with file and irrigation
- Standard luer lock syringe – one-hand operation

Pulpdent, still leading the way in calcium hydroxide preparations since developing the first pre-mixed calcium hydroxide paste in 1947.



### TempCanal Enhanced

TEK	3mL syringe + 12 endo irrigation needles (27-gauge x 1", 2-side-vent)
TE3	3mL syringe
TE4	4 x 1.2mL syringes + 20 endo irrigation needles (27-gauge x 1", 2-side-vent)

### Endodontic Irrigation Needles

27-gauge x 1" (0.4mm x 25mm), 2-side-vent

TE20N	Pkg. of 20
TE50N	Pkg. of 50
TE100N	Pkg. of 100



# PULPDENT® Endodontic Products

## PULPDENT Calcium Hydroxide Products

For Root Canal Therapy and Vital Pulp Therapy

Radiopaque, Premixed, Non-setting

Product	Root Canal Therapy	Vital Pulp Therapy	% Calcium Hydroxide	Viscosity	Contains Iodoform	Water-based	Silicone Oil-based
<b>Multi-Cal</b>	X	X	40% ± 2%	Creamy		X	
<b>TempCanal Enhanced</b>	X		40% ± 2%	Creamy		X	
<b>Pulpdent Paste</b>	X	X	40% ± 2%	Viscous		X	
<b>Forendo Paste</b>	X		22%	Creamy	X		X

MULTI	Multi-Cal	4 x 1.2mL syringes + 8 applicator tips (22-gauge x 1/2")
MULTI-3	Multi-Cal	3mL syringe
TEK	TempCanal Enhanced	3mL syringe + 12 endo irrigation needles (27-gauge x 1", 2-side-vent)
TE3	TempCanal Enhanced	3mL syringe
TE4	TempCanal Enhanced	4 x 1.2mL syringes + 20 endo irrigation needles (27-gauge x 1", 2-side-vent)
PSYK	Pulpdent Paste	3mL syringe + 24 needles (18-gauge x 1")
PSY	Pulpdent Paste	3mL syringe
FORE	Forendo Paste	2.2gr syringe + 20 applicator tips

## PULPDENT EDTA Products

Decalcifies Canal Walls, Facilitates Instrumentation, Buffered to Neutral pH

Product	Viscosity	% EDTA	Contains Lubricant	Contains Peroxide	Packaging	Recommended Delivery
<b>EDTA 17% Solution</b>	Liquid	17%			Bottle	Syringe or pipette
<b>Prep-Rite RC</b>	Viscous Paste	15%	X	X	Syringe	Pick up on file
<b>File-Rite</b>	Semi-gel	17%	X		Syringe + 30 gauge x 1" needles	Directly into canal w/ 30-ga needle

EDTA-30	EDTA 17% Solution	30mL bottle
EDTA-60	EDTA 17% Solution	60mL bottle
EDTA-120	EDTA 17% Solution	120mL bottle
EDTA480	EDTA 17% Solution	480mL bottle
PRC	Prep-Rite RC	4 x 5 gm syringes
FILE	Semi-gel	4 x 5 gm syringes + 50 needles (30-gauge x 1")

## PULPDENT Root Canal Sealer

Indicated For Permanent and Primary Teeth

For all permanent filling techniques: Pressure Syringe, Solid Core, Lentulo

Meets ANSI/ADA Specification 57 for Endodontic Filling Material

Powder: Zinc Oxide, Zinc Stearate, Calcium Phosphate, Barium Sulfate

Liquid: Eugenol, Canada Balsam

Product	Tissue Compatible	Contains Zinc and Eugenol	Radiopaque	Permanent Teeth	Primary Teeth	No Shrinkage	Slight Overfills Resorb
<b>Pulpdent Root Canal Sealer Kit</b>	X	X	X	X	X	X	X

RK	Pulpdent Root Canal Sealer Kit	15cc powder, 7.5mL liquid, mixing pad, scoop
----	--------------------------------	--

