

# HATHO THERMOPOL SET

3 Step System for easy and brilliant polishing of thermoplastic dentures.

**NEW**  
FROM  
PRESTIGE



Ref: 313203

**£52.60**



## STEP 1 PRE-POLISHING

3" HATHO multi-layer brush with 2 outer rows of Chungking bristles and 2 inner layers of Scotch-Brite™. Used with HATHO Polistar Linty-grey polishing compound made from soft wax mixture with fine grit at 1,800rpm.



## STEP 2 SHINE POLISHING

4" HATHO Cosima soft polishing buff made from 14 layers of fine, siliconized cotton cloth. Used with HATHO Polistar cream highly filled with fine grit at 1,800rpm.



## STEP 3 HIGH SHINE POLISHING

4" HATHO Mira polishing buff made from 50 layers of the finest natural cotton cloth. Used with HATHO Pro Rosa polishing compound. A handcrafted high shine polishing compound containing ultra-fine grit.

All prices exclude VAT and postage and are subject to change without prior notification.

Prestige Dental Ltd 7 Oxford Place, Bradford, West Yorkshire, BD3 0EF

Tel: 01274 721567

Email: [info@prestige-dental.co.uk](mailto:info@prestige-dental.co.uk) • [www.prestige-dental.co.uk](http://www.prestige-dental.co.uk)

If you prefer not to receive our special offers please contact us to unsubscribe



Certificate Number 8226 ISO 9001